


# GYNAECOLOGY INSTRUMENTS

A complete selection of instruments for diagnostic and surgical procedures in women's health. This category includes vaginal specula, tenacula, uterine curettes, biopsy forceps, ovum forceps, scissors, and retractors. Designed for atraumatic tissue handling, optimal access, and precise control during colposcopy, biopsy, dilatation and curettage, and major pelvic surgery. Crafted from premium stainless steel with ergonomic grips for reliable, comfortable use.



Parametrium Clamp Forceps

Hysterectomy clamp, atraumatic jaws, two gold bows



Straight

|               |            |
|---------------|------------|
| S C- G Y/07 8 | 203mm long |
| SC-GY/078a    | 233mm long |
| SC-GY/078b    | 254mm long |
| SC-GY/078c    | 304mm long |
| SC-GY/078d    | 360mm long |

Slight curve

|            |            |
|------------|------------|
| SC-GY/078e | 203mm long |
| SC-GY/078f | 233mm long |
| SC-GY/078g | 254mm long |
| SC-GY/078h | 304mm long |
| SC-GY/078i | 360mm long |

Strong curve


|            |            |
|------------|------------|
| SC-GY/078j | 203mm long |
| SC-GY/078k | 233mm long |
| SC-GY/078L | 254mm long |
| SC-GY/078n | 360mm long |

Angled

|            |            |
|------------|------------|
| SC-GY/078o | 203mm long |
| SC-GY/078p | 233mm long |
| SC-GY/078q | 254mm long |
| SC-GY/078r | 304mm long |
| SC-GY/078s | 360mm long |

Spencer-Robarts Clamp


For endometrial and pedunculated fibroid polyps, with ratchet jaw-locking mechanism, and 50° jaw opening



|           |                  |
|-----------|------------------|
| SC-GY/082 | Shaft 235mm long |
|-----------|------------------|

Museux Vulsellum Forceps

Matt finish




|           |            |
|-----------|------------|
| SC-GY/123 | 240mm long |
|-----------|------------|



Schroeder Uterine Vulsellum Forcep


1-1 teeth



|           |            |
|-----------|------------|
| SC-GY/022 | 229mm long |
|-----------|------------|

Teale Uterine Vulsellum Forcep


3-4 teeth, 230mm long



|            |          |
|------------|----------|
| SC-GY/023  | Curved   |
| SC-GY/023a | Straight |

Bonney Hysterectomy Forceps

Angled on flat, 1 - 2teeth



|           |            |
|-----------|------------|
| SC-GY/018 | 197mm long |
|-----------|------------|

Berkeley-Bonney Hysterectomy Forceps


Curve don flat, 1 - 2teeth



|           |            |
|-----------|------------|
| SC-GY/019 | 203mm long |
|-----------|------------|

Gwilliam Hysterectomy Forceps


Longitudinal serrations  
203mm long



|            |                           |
|------------|---------------------------|
| SC-GY/003a | Straight,                 |
| SC-GY/003b | Curved on flat, 1-2 teeth |

Heaney Hysterectomy Forceps


Curve don flat, 203mm long



|            |              |
|------------|--------------|
| SC-GY/016a | One tooth    |
| SC-GY/016b | Double tooth |


Maingot Hysterectomy Forceps

203mm long



|            |                             |
|------------|-----------------------------|
| SC-GY/017a | Curved on flat, 1 - 2 teeth |
| SC-GY/017b | Straight                    |

Faures Peritoneum Forceps



205mm long


|           |              |
|-----------|--------------|
| SC-GY/146 | Slight curve |
|-----------|--------------|

203mm long, longitudinal serrations

|             |                 |
|-------------|-----------------|
| SC-GY/5146  | Curved, toothed |
| SC-GY/5146a | Curved          |

Bonney Uterine Polypus Forcep

Box joint



|           |        |
|-----------|--------|
| SC-GY/020 | 1640mm |
|-----------|--------|

Noto Uterine Polypus Forceps


Matt finish



|           |        |
|-----------|--------|
| SC-GY/124 | 1640mm |
|-----------|--------|

Vant Uterine Polypus Forcep


Serrated



|           |        |
|-----------|--------|
| SC-GY/021 | 1640mm |
|-----------|--------|

Dartigues Uterine Forceps


Mattfinish



|           |        |
|-----------|--------|
| SC-GY/090 | 1670mm |
|-----------|--------|

Somer Uturus Holding Forcecep

Matt finish



|           |            |
|-----------|------------|
| SC-GY/091 | long 230mm |
| SC-GY/002 | 254mm long |

McClintock Ovum Forcecp

Curved, screwjoint



|           |            |
|-----------|------------|
| SC-GY/036 | long 241mm |
|-----------|------------|

Greenhalgh Ovum Forcecp


Screwjoint



|           |        |
|-----------|--------|
| SC-GY/037 | 1640mm |
|-----------|--------|

Keily Ovum Forceps


Matt finish



|           |        |
|-----------|--------|
| SC-GY/111 | 1620mm |
|-----------|--------|

Doyen Ovarian Forceps

Matt finish



|           |            |
|-----------|------------|
| SC-GY/118 | 210mm long |
|-----------|------------|

Wertheim Parametrium Forcep


Curved onflat, screwjoint



|           |            |
|-----------|------------|
| SC-GY/005 | 254mm long |
|-----------|------------|

Green-Armytage Haemostasis Forcep


Screw joint, 15mm jaw, 210mm long



|            |          |
|------------|----------|
| SC-GY/029a | Straight |
| SC-GY/029b | Angled   |

Kogan Endocervical Forceps


Matt finish, 275mm long



|            |         |
|------------|---------|
| SC-GY/131  | 3mm jaw |
| SC-GY/131a | 5mm jaw |
| SC-GY/131b | 7mm jaw |

Lowa Tenaculum Forceps


Matt finish, 5-6 teeth



|           |            |
|-----------|------------|
| SC-GY/094 | 260mm long |
|-----------|------------|

Skene Tenaculum Forceps


Bright finish



|           |            |
|-----------|------------|
| SC-GY/120 | 240mm long |
|-----------|------------|

Pratt Tenaculum Forceps

Matt finish



|           |            |
|-----------|------------|
| SC-GY/095 | 280mm long |
|-----------|------------|

Coil Removing Forcep


With serrated rounded jaws  
4mm wide



|           |                   |
|-----------|-------------------|
| SC-GY/075 | 160mm to shoulder |
|-----------|-------------------|

Willet Martel Forceps

Matt finish




|           |            |
|-----------|------------|
| SC-GY/121 | 190mm long |
|-----------|------------|



Littlewood Tissue Forceps


2-3 teeth, matt finish



|           |            |
|-----------|------------|
| SC-GY/275 | 190mm long |
|-----------|------------|

Eppendorf Biopsy Punch Forcep

200mm to shoulder




Bite shown actual size

|                  |                |
|------------------|----------------|
| SC-GY/14-314-200 | 3mm x 8mm bite |
|------------------|----------------|

Kevorkian Biopsy Punch Forcep

255mm to shoulder




Bite shown actual size

|                  |                |
|------------------|----------------|
| SC-GY/14-329-250 | 3mm x 8mm bite |
|------------------|----------------|

Tischler Morgan Biopsy Forcep

210mm to shoulder



Bite shown actual size

|                  |                |
|------------------|----------------|
| SC-GY/14-328-220 | 3mm x 7mm bite |
|------------------|----------------|

Tischler Morgan Biopsy Punch Forcep

295mm to shoulder




Bite shown actual size

|           |                |
|-----------|----------------|
| SC-GY/088 | 2mm x 7mm bite |
|-----------|----------------|

Baby Tischler Biopsy Forcep

210mm to shoulder



 Bite shown actual size

|                  |                |
|------------------|----------------|
| SC-GY/14-325-250 | 2mm x 5mm bite |
|------------------|----------------|

Baby Tischler Biopsy Punch Forcep

295mm to shoulder




 Bite shown actual size

|           |                |
|-----------|----------------|
| SC-GY/087 | 2mm x 5mm bite |
|-----------|----------------|

Van Doren Biopsy Forceps

Bright finish, 260mm overall length




 Bite shown actual size

|           |                |
|-----------|----------------|
| SC-GY/112 | 5mm x 9mm bite |
|-----------|----------------|



Van Doren Cervical Biopsy Punch Forcep

With curved basket jaw, angled shanks



|           |            |
|-----------|------------|
| SC-GY/025 | 240mm long |
|-----------|------------|

Wittner Cervical Biopsy Punch Forcep


240mm longshaft



|            |                       |
|------------|-----------------------|
| SC-GY/024a | Jaws curved to side   |
| SC-GY/024b | Jaws straight to side |

Barnes Obstetric Forcep


With Simpson handles and Neville axis-traction rod



|           |            |
|-----------|------------|
| SC-GY/027 | 381mm long |
|-----------|------------|

Haig Ferguson Obstetric Forcep


Tape slots in blades and axis traction rod



|           |            |
|-----------|------------|
| SC-GY/028 | 356mm long |
|-----------|------------|

Kieland Obstetric Forcep


Matt finish



|           |            |
|-----------|------------|
| SC-GY/089 | 410mm long |
|-----------|------------|

Simpson Braun Obstetric Forcep


Matt finish



|           |            |
|-----------|------------|
| SC-GY/086 | 380mm long |
|-----------|------------|

Wrigley Obstetric Forcep


Short model for caesarian section



|           |            |
|-----------|------------|
| SC-GY/026 | 279mm long |
|-----------|------------|

Auvard Vaginal Speculum

Weighted



|                |                       |
|----------------|-----------------------|
| S C- G Y/O 07a | Small blade 38 x 90mm |
| SC-GY/007b     | Large blade 51 x 90mm |

Graves Vaginal Speculum


Bright finish



|           |            |
|-----------|------------|
| SC-GY/103 | 115 x 35mm |
|-----------|------------|


Cusco Vaginal Speculum

Stainless steel



| Insulated  |                           | With Suction Tube |                            | Black Finish |                            |
|------------|---------------------------|-------------------|----------------------------|--------------|----------------------------|
| SC-GY/015d | Extra small - 98mm x 16mm | SC-GY/134         | Small - 105mm x 25mm       | SC-GY/015k   | Small - 105mm x 25mm       |
| SC-GY/015c | Small - 105mm x 25mm      | SC-GY/135         | Medium - 107mm x 29mm      | SC-GY/015j   | Medium - 107mm x 29mm      |
| SC-GY/015f | Medium - 107mm x 29mm     | SC-GY/136         | Large - 115mm x 35mm       | SC-GY/015i   | Large - 115mm x 35mm       |
| SC-GY/015e | Large - 115mm x 35mm      | SC-GY/137         | Extra large - 135mm x 38mm | SC-GY/015h   | Extra large - 135mm x 38mm |


Collin Vaginal Speculum



|            |                              |
|------------|------------------------------|
| SC-GY/014a | Small - blades 96mm x 25mm   |
| SC-GY/014b | Medium - blades 108mm x 30mm |
| SC-GY/014c | Large - blades 108mm x 35mm  |

Nott Speculum


Stainless steel



|             |                             |
|-------------|-----------------------------|
| SC-GY/04102 | Blade 110mm long, 30mm wide |
|-------------|-----------------------------|

Pederson Vaginal Speculum


Stainless steel



|               |                   |
|---------------|-------------------|
| S C- G Y/07 B | 120mm x 25mm wide |
| SC-GY/078a    | 90mm x 20mm wide  |

Winterton Vaginal Speculum


Stainless steel



|              |                             |
|--------------|-----------------------------|
| S C- G Y/074 | Blade 115mm long, 25mm wide |
|--------------|-----------------------------|

Modified Speculum


Stainless steel, for elderly patients



|                |                             |
|----------------|-----------------------------|
| S C- G Y/074 a | Blade 180mm long, 18mm wide |
|----------------|-----------------------------|

Sims Duckbill Speculum


Double ended



|            |                                   |
|------------|-----------------------------------|
| SC-GY/012a | Small - 27mm x 29mm x 130mm long  |
| SC-GY/012b | Medium - 30mm x 33mm x 158mm long |
| SC-GY/012c | Large - 33mm x 38mm x 184mm long  |

Sims Duckbill Speculum


Double ended, insulated



|            |                                   |
|------------|-----------------------------------|
| SC-DIA/031 | Small - 27mm x 29mm x 130mm long  |
| SC-DIA/032 | Medium - 30mm x 33mm x 158mm long |
| SC-DIA/033 | Large - 33mm x 38mm x 184mm long  |


Breisky Vaginal Retractor

Bright finish, matt handle  
320mm long




|           |            |
|-----------|------------|
| SC-GY/110 | 25mm blade |
|-----------|------------|

Doyen Rectangular Retractor




|            |                    |
|------------|--------------------|
| SC-GY/008a | Small - 48 x 60mm  |
| SC-GY/008b | Medium - 48 x 89mm |
| SC-GY/008c | Large - 48 x 124mm |

Doyen Retractor




|            |                   |
|------------|-------------------|
| SC-GY/009a | Blade 95mm x 27mm |
|------------|-------------------|

Doyen Retractor



|            |                    |
|------------|--------------------|
| SC-GY/009b | Blade 100mm x 30mm |
|------------|--------------------|

St Mark's Hospital Pattern Retractor




With lip

|           |                    |
|-----------|--------------------|
| SC-GY/065 | Blade 125mm x 60mm |
| SC-GY/066 | Blade 175mm x 60mm |

Without lip

|           |                    |
|-----------|--------------------|
| SC-GY/067 | Blade 125mm x 60mm |
| SC-GY/068 | Blade 175mm x 60mm |

Wertheim Retractor



|           |                   |
|-----------|-------------------|
| SC-GY/010 | Blade 95mm x 38mm |
|-----------|-------------------|


Landon Vaginal Retractor



| SC-GY/011 | Blade 89mm x 25mm |
|-----------|-------------------|
|-----------|-------------------|

Lateral Vaginal Retractor

Blades bright finish, matt handle



| SC-GY/128  | 60 x 20mm blade  |
|------------|------------------|
| SC-GY/128a | 75 x 22mm blade  |
| SC-GY/128b | 110 x 23mm blade |

AVAILABLE INSULATED  
ADD SUFFIX 'INS'

Whitlow Retractor

Vaginal retractor, for use with supplied blades. Points of retraction move around almost 360°. Supplied complete with sterilising basket and three sets of blades.



Frame

225mmat widest point

Blades

1 x anterior and 1 x posterior 90°, 25mm x 100mm  
2 x lateral 90°, 25mm x 70mm  
2 x lateral 75°, 25mm x 70mm


| SC-GY/084 |
|-----------|
|-----------|



Kevorkian Biopsy Curette

313mm long

| SC-GY/105 | With basket    |
|-----------|----------------|
| SC-GY/127 | Without basket |



Novak Curette

| SC-GY/119 | 230mm long |
|-----------|------------|
|-----------|------------|



Recamier Uterine Curette

Matt finish, 310mm long

| SC-GY/108 | Blunt |
|-----------|-------|
| SC-GY/109 | Sharp |



Rheinstadter Flushing Curette

290mm long

Sharp

| SC-GY/076a | Sharp, 8mm  |
|------------|-------------|
| SC-GY/076b | Sharp, 14mm |

Blunt

| S C- G Y/076 c | Blunt, 8mm  |
|----------------|-------------|
| SC-GY/076d     | Blunt, 14mm |



Sims Uterine Curette


Single ended, 255mm long, matt finish

Sharp

| SC-GY/075a | Small 6mm   |
|------------|-------------|
| SC-GY/075b | Medium 10mm |
| SC-GY/075c | Large 13mm  |

Blunt

| SC-GY/075d | Small 6mm   |
|------------|-------------|
| SC-GY/075e | Medium 10mm |
| SC-GY/075f | Large 13mm  |



Sims Uterine Curette

Double ended, 250mm long

Sharp

| SC-GY/044 | Small 7 x 5mm   |
|-----------|-----------------|
| SC-GY/045 | Medium 10 x 7mm |
| SC-GY/046 | Large 13 x 10mm |


Blunt

| S C- G Y/0 47 | Small 7 x 5mm   |
|---------------|-----------------|
| SC-GY/048     | Medium 10 x 7mm |
| SC-GY/049     | Large 13 x 10mm |



Rotunda Hospital Pattern Catheter Catheter

With solid end, chromium plated, 203mm long      Oval, solid end, chromium plated



|            |         |
|------------|---------|
| SC-GY/035a | 8 E.G.  |
| SC-GY/035b | 10 E.G. |



|           |            |
|-----------|------------|
| SC-GY/034 | 165mm long |
|-----------|------------|

Catheter Introducer

457mm long



|               |                 |
|---------------|-----------------|
| S C- G Y/07 1 | Straight, 2mm ø |
| SC-GY/071a    | Curved, 2mm ø   |

Catheter Introducer



|            |                  |
|------------|------------------|
| SC-GY/071b | Fine, 395mm long |
|------------|------------------|

Catheter Introducer

Bulbous end, with eye in handle



|            |                          |
|------------|--------------------------|
| SC-GY/071c | Curved, 2mm x 410mm long |
|------------|--------------------------|

Catheter Introducer

With eye in handle, curved



|            |                  |
|------------|------------------|
| SC-GY/071d | 2mm x 420mm long |
|------------|------------------|

Catheter Introducer

Fixed handle



|            |            |
|------------|------------|
| SC-GY/071e | 445mm long |
|------------|------------|

Catheter Introducer



|            |            |
|------------|------------|
| SC-GY/071f | 400mm long |
|------------|------------|


Catheter Introducer

410mm long



|            |            |
|------------|------------|
| SC-GY/071g | 410mm long |
|------------|------------|


Barnes Episiotomy Scissors



|           |            |
|-----------|------------|
| SC-GY/032 | 127mm long |
|-----------|------------|

Braun Stadler Episiotomy Scissors

Matt finish



|        |       |
|--------|-------|
| SC-    | 140mm |
| GY/063 | long  |
| SC-    | 222mm |
| GY/064 | long  |

Vant Episiotomy Scissors



|           |                         |
|-----------|-------------------------|
| SC-GY/033 | 76mm blades, 241mm long |
|-----------|-------------------------|

Busch Umbilical Scissors

Matt finish



|           |            |
|-----------|------------|
| SC-GY/116 | 160mm long |
|-----------|------------|

Umbilical Scissors


American pattern



|           |            |
|-----------|------------|
| SC-GY/031 | 102mm long |
|-----------|------------|

Simms Uterine Scissors


Curved,tungstencarbide,withtwogold bows



|           |            |
|-----------|------------|
| SC-GY/115 | 200mm long |
|-----------|------------|

Dubois Uterine Scissors

Matt finish



|           |            |
|-----------|------------|
| SC-GY/114 | 270mm long |
|-----------|------------|

Uterine Scissors


Curved on flat



|            |            |
|------------|------------|
| SC-GY/030a | 203mm long |
|------------|------------|

Uterine Scissors



Angled on flat



|            |            |
|------------|------------|
| SC-GY/030b | 230mm long |
|------------|------------|

Parametrium Scissors

225mmlong



Tungsten carbide, standard, with two gold bows

|            |              |
|------------|--------------|
| SC-GY/080  | Slight curve |
| SC-GY/080a | Strong curve |
| SC-GY/080b | Acute curve  |


Tungsten carbide, supercut, with one gold and one black bow

|            |              |
|------------|--------------|
| SC-GY/081  | Slight curve |
| SC-GY/081a | Strong curve |
| SC-GY/081b | Acute curve  |

Miscellaneous

BonneyMyomectomy Clamp

With 51mm jaws, screw joint



|           |            |
|-----------|------------|
| SC-GY/004 | 267mm long |
|-----------|------------|

Sims Vaginal Depressor

Doupleended



|           |            |
|-----------|------------|
| SC-GY/039 | 260mm long |
|-----------|------------|

Spackman Uterine Cannula

Chromeplated,360mmlong



|            |                          |
|------------|--------------------------|
| SC-GY/040  | 12 F.G. with rubber cone |
| SC-GY/040a | 14 F.G. with metal cone  |
| SC-GY/040b | 14 F.G. with rubber cone |
| SC-GY/040c | 18 F.G. with metal cone  |
| SC-GY/040d | 18 F.G. with rubber cone |

Doyen Myoma Screw Large



|            |            |
|------------|------------|
| SC-GY/006a | 178mm long |
|------------|------------|


Doyen Myoma Screw Small



|            |            |
|------------|------------|
| SC-GY/006b | 178mm long |
|------------|------------|


Lake Seeker

Double ended, for vestibule



|            |            |
|------------|------------|
| SC-ENT/040 | 184mm long |
|------------|------------|

Lempert Rugine



|             |                    |
|-------------|--------------------|
| SC-ENT/064a | Narrow, 156mm long |
| SC-ENT/065a | Broad, 171mm long  |

Irwin-Moore Tonsil Needle


203mmlong



|             |                 |
|-------------|-----------------|
| SC-ENT/164a | Curved to right |
| SC-ENT/164a | Curved to left  |


Negus Pusher

“Knottier”, metal handle



|            |            |
|------------|------------|
| SC-ENT/165 | 203mm long |
|------------|------------|


Stanford Cade Insufflator



|            |            |
|------------|------------|
| SC-ENT/188 | 175mm long |
|------------|------------|

Tubing Mount

Foruse with catheters, Luer



|            |      |
|------------|------|
| SC-ENT/016 | 25mm |
|------------|------|

Galabin Uterine Sound

Calibrated 4-24, bright finish



|           |            |
|-----------|------------|
| SC-GY/106 | 305mm long |
|-----------|------------|

Sims Uterine Sound

Calibrated 4-24, bright finish



|           |            |
|-----------|------------|
| SC-GY/107 | 320mm long |
|-----------|------------|

Umbilical Dilator



|           |           |
|-----------|-----------|
| SC-GY/038 | 95mm long |
|-----------|-----------|

Hawkin Ambler Cervix Dilators

Conical, 171mm long



|           |               |
|-----------|---------------|
| SC-GY/050 | Set of 12     |
| SC-GY/051 | Size 3 - 6mm  |
| SC-GY/052 | Size 4 - 7mm  |
| SC-GY/053 | Size 5 - 8mm  |
| SC-GY/054 | Size 6 - 9mm  |
| SC-GY/055 | Size 7 - 10mm |
| SC-GY/056 | Size 8 - 11mm |

|           |                |
|-----------|----------------|
| SC-GY/057 | Size 9 - 12mm  |
| SC-GY/058 | Size 10 - 13mm |
| SC-GY/059 | Size 11 - 14mm |
| SC-GY/060 | Size 12 - 15mm |
| SC-GY/061 | Size 13 - 16mm |
| SC-GY/062 | Size 14 - 17mm |

Hegar Cervix Dilator

Single ended, 185mm long




|             |       |
|-------------|-------|
| SC-GY/70SE1 | 1mm ø |
| SC-GY/70SE2 | 2mm ø |
| SC-GY/70SE3 | 3mm ø |
| SC-GY/70SE4 | 4mm ø |
| SC-GY/70SE5 | 5mm ø |
| SC-GY/70SE6 | 6mm ø |
| SC-GY/70SE7 | 7mm ø |
| SC-GY/70SE8 | 8mm ø |
| SC-GY/70SE9 | 9mm ø |

|              |        |
|--------------|--------|
| SC-GY/70SE10 | 10mm ø |
| SC-GY/70SE11 | 11mm ø |
| SC-GY/70SE12 | 12mm ø |
| SC-GY/70SE13 | 13mm ø |
| SC-GY/70SE14 | 14mm ø |
| SC-GY/70SE15 | 15mm ø |
| SC-GY/70SE16 | 16mm ø |
| SC-GY/70SE17 | 17mm ø |

Hegar Cervix Dilator

Double ended, 200mm long



|                |                                               |
|----------------|-----------------------------------------------|
| S C- G Y/070 j | Set of 8: from sizes 1 and 2mm to 15 and 16mm |
| SC-GY/070h     | 1 and 2mm ø                                   |
| S C- G Y/070   | 3 and 4mm ø                                   |
| SC-GY/070a     | 5 and 6mm ø                                   |
| S C- G Y/070 b | 7 and 8mm ø                                   |

|                |               |
|----------------|---------------|
| S C- G Y/070 c | 9 and 10mm ø  |
| SC-GY/070d     | 11 and 12mm ø |
| SC-GY/070e     | 13 and 14mm ø |
| SC-GY/070f     | 15 and 16mm ø |
| SC-GY/070g     | 17 and 18mm ø |