


OPHTHALMIC/ MICRO

A delicate, high-precision collection for ocular and microsurgical procedures. Includes micro forceps, scissors, needle holders, cannulas, speculums, and knives designed for fine tissue manipulation under high magnification. Ultra-fine tips, smooth spring-action handles, and non-glare surfaces provide exceptional control and visibility in confined fields. Ideal for anterior/posterior segment surgery, corneal repair, and delicate micro-anastomosis. Crafted from premium stainless steel for lasting sharpness and tactile sensitivity.



Iris Forceps


1-2teeth, 114mm long



SC-OP/01	Straight
SC-S-C-O-P/02	Curved

Iris Forceps


Serrated, 114mm long



SC-OP/03	Straight
SC-OP/04	Curved

Iris Forceps


Extra fine, serrated, 94mm long



SC-OP/06	Curved
----------	--------

Iris Fixation Forceps

Extra fine, 1-2 teeth, 102mm long



SC-OP/05	Curved
----------	--------

Barraquer Corneal Forceps


76mm long



SC-OP/023	1-2 teeth
-----------	-----------

Epilation Forceps


89mm long



SC-OP/07	Oblique jaw
----------	-------------

Epilation Forceps


89mm long



SC-OP/08	Straight ends, 3mm wide
SC-OP/08a	Round ends, 3mm wide

Suture/Dissecting Forceps


102mm long



SC-OP/09	Extra fine
----------	------------

Barraquer Fixation Forceps


102mm long



SC-OP/09a	9-10 teeth
-----------	------------

Jameson Evans Cilia Forceps


120mm long



SC-OP/012	
-----------	--

Tarsal Cyst Forceps


89mm long



SC-OP/015	Round ends
-----------	------------

St Martin Corneal Suturing Forceps


1-2 teeth, 89mm long



SC-OP/022	For fine work
-----------	---------------

St Martin Micro Corneal Suturing Forceps


1-2 teeth, 89mm long



SC-OP/022a	Fine pattern
------------	--------------

Castroviejo Corneal Suture Forceps


94mm long



SC-OP/013	1-2 teeth
-----------	-----------

Arruga Intra Capsule Forceps


For removal of residual lens material during extra-capsula extraction, 120mm long



SC-OP/024	Curved tips
-----------	-------------

Desmarre Entropian Forceps


89mm long



SC-OP/014	Right eye
SC-OP/014a	Left eye


Wilde Entropian Forceps

3sizes, 89mm long



SC-OP/015a	With screw action
------------	-------------------


de Wecker Iris Scissors



SC-OP/031	114mm long
-----------	------------

Castroviejo Scissors

102mm long



SC-OP/025	Curved
-----------	--------


Barraquer Iris Scissors



SC-OP/031a	62mm long
------------	-----------

Iris Scissors


102mm long



SC-OP/032	Straight
SC-OP/033	Curved

Pooley Conjunctival Scissors


89mm long



SC-OP/017	For fine dissection
-----------	---------------------

Thomas Strabismus Scissors

114mm long



SC-OP/029	Straight
SC-OP/030	Curved

Strabismus Scissors


Super cut, 114mm long



SC-GI/475	Straight
SC-GI/475a	Curved

Iris Scissors

Super cut, 114mm long



SC-O P/O 476	Straight
SC-O P/O 476 a	Curved

Westcott Scissors

120mm long

SC-OP/251	Curved, blunt
SC-OP/249	Curved, sharp
SC-OP/250	Straight, blunt
SC-OP/248	Straight, sharp


Vannas Scissors

Scissors

SC-OP/241	Blunt, curved 45°, 110mm
SC-OP/245	Curved 85mm
SC-OP/244	Straight 85mm
SC-OP/246	Strongly curved 85mm


Micro Scissors

SC-OP/278	Delicate, straight, 85mm
SC-OP/280	Delicate, curved, 85mm
SC-OP/279	Delicate, strongly curved in plane, 85mm



Castroviejo Needle Holder

SC-OP/021	140mm long, with catch
SC-OP/020	140mm long, without catch



Diamond dust jaws


SC-THO/060	180mm long, without catch
SC-THO/061	200mm long, without catch
SC-THO/062	215mm long, without catch



ALSO AVAILABLE WITH CATCH - ADD SUFFIX CTH TO

Tungsten carbide jaws

SC-THO/058a	180mm long
SC-THO/058	200mm long
SC-THO/059	215mm long




ALSO AVAILABLE WITH CATCH - ADD SUFFIX CTH TO

Silcock Needle Holder

140mm long


SC-OP/21a	With catch
SC-OP/20	Without catch



Micro Needle Holder

Round knurled handle, 120mm long


SC-OP/YA4020	Curved, without catch
SC-OP/YA4050	Curved, with catch
SC-OP/WA2010	Straight, without catch



Micro Needle Holder


Curved, without catch

SC-OP/YA4010	110mm long
--------------	------------



Lang Eyelid Retractor

76mm long



SC-OP/011b	Left eye
SC-OP/011c	Right eye

Desmarre Eyelid Retractor


Single, 140mm long



SC-OP/013	12.7mm wide blade
SC-OP/014	15.9mm wide blade
SC-OP/015	19mm wide blade

Barraquer Eye Speculum


Stainlesssteelwire,44mm long



SC-OP/006	Standard
SC-OP/007	Small

Speculum

Wire, 30mm wide x 45mm long



SC-OP/8004012	Titanium
SC-OP/8004012ss	Stainless

Speculum

Wire,25mmwidex 40mm long



SC-OP/8004013	Titanium
SC-OP/8004013ss	Stainless

Speculum


Wire,20mmwidex 30mm long



SC-OP/8004013a	Titanium
SC-OP/8004013a-ss	Stainless

Clark Speculum

Wire, 35mm wide x 84mm long



SC-SH/5746L	Infant - left
SC-SH/5746R	Infant - right

Lang Lachrymal Cannula

38mm long



SC-OP/026	Cannula = 1mm ø
-----------	-----------------

Cannula

Curved,44mm long



SC-OP/026	Flat ended
-----------	------------

Harrison Butler Lachrymal Cannula

38mm long



SC-OP/027	Cannula = 0.6mm ø
-----------	-------------------