




NEURO/ SPINAL SURGERY

A precision-focused collection for delicate cranial and spinal procedures. Includes self-retaining retractors, rongeurs (Kerrison, pituitary), dissectors, microsurgical forceps, scissors, elevators, and curettes. Designed for minimal tissue trauma, precise bone removal, and controlled access to confined neural structures. Fine-tipped, balanced instruments with non-glare finishes support accurate manipulation under the microscope or loupes. Crafted from durable, high-grade stainless steel for reliable long-term performance.



Mastoid Retractor


3 x 3 blunt hooks



SC-NEU/N30350	114mm long
---------------	------------

Jansen (Shea) Retractor


Self retaining, blunt blades



SC-NEU/028	114mm long
------------	------------

Cone Retractor


Selfretaining,hingedarms, for use in ventriculography



SC-NEU/027	165mm long
------------	------------

Jefferson Retractor


Self retaining, curved, blunt pointed blades



SC-NEU/029	140mm long
------------	------------

Henderson Retractor


Self retaining, straight, 4x4 blunt pointed blades



SC-NEU/030a	178mm long
-------------	------------

Cairns Scalp Retractor


Curved, self retaining, swivel blades



SC-NEU/025b	140mm long
-------------	------------

Schnitker Scalp Retractor


Curved, self retaining, blunt swivel blades



SC-NEU/026	140mm long
------------	------------

Cone Laminectomy Retractor

Hinged arms, with cam rack




SC-NEU/051	254mm long
------------	------------

Cone Laminectomy Retractor

Hinged arms


SC-NEU/052a	Left, 165mm long
SC-NEU/052b	Right, 165mm long
SC-NEU/052c	Left, 254mm long
SC-NEU/052d	Right, 254mm long



Harvey-Jackson Laminectomy Retractor

Hinged arms, blades 38mm wide x 44mm deep


SC-NEU/053	305mm long
------------	------------



Mayo-Adson (Beckmann) Laminectomy Retractor

Hinged arms

SC-NEU/054	292mm long
------------	------------




Northfield (Cairns) Retractor

Hinged arms, one pair of blades

4 prong : 25mm wide x 16mm deep
and one pair of blades


5 prong: 32mm wide x 25mm deep

SC-NEU/055	292mm long
------------	------------




Meyerding Hemi-Laminectomy Retractor

SC-NEU/216	18cm long 7" blade size 50 x 10mm
SC-NEU/219	18cm long 7" blade size 50 x 20mm
SC-NEU/217	18cm long 7" blade size 60 x 10mm
SC-NEU/220	18cm long 7" blade size 60 x 20mm
SC-NEU/218	18cm long 7" blade size 70 x 10mm
SC-NEU/221	18cm long 7" blade size 70 x 20mm
SC-NEU/224	18cm long blade left 90 x 10mm
SC-NEU/225	18cm long blades left 90 x 20mm
SC-NEU/222	18cm long blades right 90 x 10mm
SC-NEU/223	18cm long blades right 90 x 20mm



Cairns-Northfield Laminectomy Retractor

203mm arm length with three pairs of interchangeable swivelling blades, 76mm, 64mm, and 51mm



SC-NEU/036

Cairns 'Rake' Scalp Retractor

Sharp, 210mm long



SC-NEU/032a	3 prong, 16mm wide
SC-NEU/032b	4 prong, 19mm wide

Cairns (Cushing) Scalp Retractor


Solid blade



SC-NEU/033	16mm wide, 210mm long
------------	-----------------------

Cairns Strip Retractor

Malleable, 203mm long
8 sizes: 6, 8, 10, 11, 13, 16, 19, 25mm



SC-NEU/034a	20 SWG
SC-NEU/034b	18 SWG

Hibbs Retractor

203mm long



SC-NEU/050a	Blades 10mm x 38mm
SC-NEU/050b	Blades 16mm x 54mm

Krayenbuhl Retractor

6mm wide



SC-NEU/5905	210mm long
SC-NEU/5905a	240mm long

Love Retractor

6mm wide



SC-NEU/5905b	210mm long
--------------	------------

Love Retractor

203mm long



S C- N E U/07 7a	Straight, blade 5mm x 5mm
SC-NEU/077b	Shaft angled 90° at handle
S C- N E U/07 7c	Shaft angled 45° at handle

O'Connell Retractor


Double ended



S C- N E U/076	216mm long
----------------	------------

Cobb Spinal Elevator

13mm wide, knurled handle



Cobb Spinal Elevator

With blue acetal handles, 279mm long

SC-NEU/073a	10mm wide
SC-NEU/073b	13mm wide
SC-NEU/073c	19mm wide

ALSO AVAILABLE WITH ACETAL SUREGRIP HANDLE
ADD SUFFIX BH TO CATALOGUE NUMBER



Adson Periosteal Elevator

Serrated grip, 165mm long



straight

SC-OR/N028	7mm wide
------------	----------

curved

SC-OR/N029	5mm wide
SC-OR/N030	7mm wide

Cushing Periosteal Elevator

Slightly curved end



SC-NEU/071	197mm long
------------	------------

McKissock Spade Periosteal Elevator


Straight



SC-NEU/072	197mm long
------------	------------

Pennybacker Elevator

Curved, sharp, 171mm long



SC-OR/N025	6mm wide
------------	----------

Cairns Guide and Director

Malleable



SC-NEU/080	127mm long
------------	------------

Dissector

200mm long



SC-NEU/N30210	Soft
SC-NEU/N30220	Soft, with hole

Penfield Dissector



No 1

Broad curve, with 7.5mm sharp round spoon, 184mm long




No 2

7.5mm sharp round spoon with wax packer, 178mm long



No 3

Fullcurve and wax packer, 197mm long



No 4

Slightcurve, 4mm wide, 216mm long

SC-NEU/069a	No 1
SC-NEU/069b	No 2
SC-NEU/069c	No 3
SC-NEU/069d	No 4

Pennybacker Probe Dissector


Double ended



SC-NEU/101	229mm long
------------	------------

Swedish Dissector

Double ended



SC-NEU/100a	Light model, 178mm long
SC-NEU/100b	Heavy model, 200mm long

Watson Cheyne Dissector

Fine, double-ended, with probe end



SC-GI/287	127mm long
SC-GI/288	191mm long
SC-GI/288a	230mm long
SC-GI/288b	255mm long

Hooks

Adson Hook


178mm long



SC-NEU/078a	Angled
SC-NEU/078b	Straight
SC-NEU/078b-b	Straight, blunt

Adson Oliver Neuro Hook

191mm long



SC-NEU/082a	Angled downward, blunt
SC-NEU/082b	Sharp
SC-NEU/082c	Angled upward, blunt

Cairns Tissue Hook


127mm long



SC-NEU/081b	Half round, sharp
SC-NEU/081c	Angled, sharp

Oliver Tissue Hook

Blunt




SC-NEU/081a	127mm long
-------------	------------

Curettes

Falconer Scoop Curette

Double ended ,angled ends



SC-NEU/099	191mm long
------------	------------

Adson Frazier Suction Tube


Angled, fully insulated, 180mm long



SC-NEU/094c	8 F.G.
-------------	--------

Retractor Suction Tube


With stilette, 220mm long



SC-NEU/SK001STIL	Replacement stilette
------------------	----------------------

Frazier Suction Tubes


Malleable, silver, angled, 180mm long



SC-NEU/092a	7 F.G.
SC-NEU/092b	8 F.G.
SC-NEU/092c	9 F.G.

Dandy Rongeur


Straight



SC-NEU/134	3mm bite, 178mm shaft
------------	-----------------------

Falconer Rongeur


Oblique jaws



SC-NEU/122a	Upward cutting, 3mm bite, 110mm shaft
SC-NEU/122b	Upward cutting, 5mm bite, 110mm shaft
SC-NEU/122c	Downward cutting, 3mm bite, 85mm shaft
SC-NEU/122d	Downward cutting, 5mm bite, 110mm shaft

Beck Rongeur

Straight, for intervertebral disc



SC-NEU/125	2mm bite, 130mm shaft
------------	-----------------------

Cushing Rongeur


130mm shaft



SC-NEU/127-2	2mm bite
SC-NEU/127-3	3mm bite
SC-NEU/127-4	4mm bite
SC-NEU/127-5	5mm bite
SC-NEU/127-6	6mm bite

Cushing Rongeur

130mm shaft



150° upward cutting


SC-NEU/128c	2mm bite
SC-NEU/128	3mm bite
SC-NEU/128a	4mm bite
SC-NEU/128b	5mm bite

150° downward cutting

SC-NEU/130c	2mm bite
SC-NEU/130	3mm bite
SC-NEU/130a	4mm bite
SC-NEU/130b	5mm bite

Love-Grunwald Rongeur

Cup shaped, jaws 3mm x 10mm, 178mm shaft



SC-NEU/131a	Straight
SC-NEU/131b	Curved upward
SC-NEU/131c	Curved downward

Mackay Rongeur


140°, straight, 114mm shaft



S C-S H /424 9	5.5mm
----------------	-------

Northfield Rongeur

Curved, upward cutting



SC-NEU/132	4mm bite, 133mm shaft
------------	-----------------------

Pennybacker Rongeur

120mm shaft



SC-NEU/133	3mm bite
SC-NEU/133-4	4mm bite

St Bart's Rongeur


114mm shaft



SC-ENT/140a	3mm bite
-------------	----------

St Bart's Rongeur

140°, downward cutting, 114mm shaft



SC-ENT/140ac	3mm
--------------	-----

St Bart's Rongeur


140°, upward cutting, 114mm shaft



SC-ENT/140b	3mm
-------------	-----

Daniel Rongeur

Miniature, for neuro pathology



SC-NEU/104	1.5mm bite, 127mm long
------------	------------------------

Cairns Rongeur


Finejaws, angled on flat, serrated jaws, curved handles



SC-NEU/103	4mm bite 152mm long
------------	---------------------

Danish Pattern Rongeur


Angled on flat, fine jaws, curved handles, screw joint



SC-NEU/109	4mm bite, 203mm long
------------	----------------------

Falconer Rongeur


Angled to side, light, screw joint



SC-NEU/110	6mm bite, 203mm long
------------	----------------------

Wilms Rongeur


Angled on flat, narrow jaws



SC-NEU/114	4mm bite, 203mm long
------------	----------------------

Glasgow Pattern Rongeur


Slightly curved on flat, screw joint



SC-NEU/111	4mm bite, 146mm long
------------	----------------------

Heanley Rongeur


Long finest raight jaws, screw joint



SC-NEU/112	4mm bite, 178mm long
------------	----------------------

Jansen Rongeur


Shortcurve, screwjoint



SC-NEU/105	8mm bite, 152mm long
------------	----------------------

Pennybacker Rongeur


Straight, fine jaws, 203mm long



SC-NEU/107	4mm bite
------------	----------

Northfield Rongeur


Straight, 10mm bite



SC-NEU/106	One shallow blade, 200mm long
------------	-------------------------------

Stookey Rongeur


Angled to side, fine jaws



SC-NEU/108	3mm bite, 178mm long
------------	----------------------

Trotter Rongeur


Angled to side, strong, screw joint




SC-NEU/113	10mm bite, 203mm long
------------	-----------------------

Luer-Jansen Rongeur

Light, curved on flat, compound action



SC-NEU/116	4mm bite, 197mm long
------------	----------------------




Olivecrona Rongeur

Angled to side, compound action, 229mm long



SC-NEU/119a	5mm bite
SC-NEU/119b	6mm bite




Ruskin Rongeur

Curved on flat, compoundaction



SC-NEU/118	6mm bite, 184mm long
------------	----------------------




Leksell Rongeur

Curved to side, compound action



SC-NEU/115a	Small, 6mm bite, 229mm long
SC-NEU/115b	Large, 8mm bite, 254mm long




Luer Stille Bone Rongeur

Curved, 250mm long




SC-OR/071b	11mm bite
SC-NEU/117	Heavy, 229mm long



Swedish Pattern Laminectomy Shears

Double crank, compound action

SC-NEU/121	267mm long
------------	------------



Rongeurs


Kerrison Rongeurs

Clowarddisc punch (notillustrated)
130° upward cutting, 190mm long
BLACK CERAMIC COATED KERRISONS AVAILABLE
PLEASE ADD SUFFIX -BK WHEN ORDERING

SC-NEU/CL002	2mm wide
SC-NEU/CL003	3mm wide
SC-NEU/CL004	4mm wide
SC-NEU/CL005	5mm wide
SC-NEU/CL006	6mm wide



Kerrison punch



Non-Detachable

90°upwardcutting

SC-NEU/136	180mm(L) x 1mm(W)
SC-NEU/137	180mm(L) x 2mm(W)
SC-NEU/138	180mm(L) x 3mm(W)
SC-NEU/139	180mm(L) x 4mm(W)
SC-NEU/140	180mm(L) x 5mm(W)
SC-NEU/141	230mm(L) x 2mm(W)
SC-NEU/142	230mm(L) x 3mm(W)
SC-NEU/143	230mm(L) x 4mm(W)
SC-NEU/144	230mm(L) x 5mm(W)

90° downward cutting

SC-NEU/154	180mm(L) x 1mm(W)
SC-NEU/155	180mm(L) x 2mm(W)
SC-NEU/156	180mm(L) x 3mm(W)
SC-NEU/157	180mm(L) x 4mm(W)
SC-NEU/158	180mm(L) x 5mm(W)
SC-NEU/159	230mm(L) x 2mm(W)
SC-NEU/160	230mm(L) x 3mm(W)
SC-NEU/161	230mm(L) x 4mm(W)
SC-NEU/162	230mm(L) x 5mm(W)

130° upward cutting

SC-NEU/172	180mm(L) x 1mm(W)
SC-NEU/173	180mm(L) x 2mm(W)
SC-NEU/174	180mm(L) x 3mm(W)
SC-NEU/175	180mm(L) x 4mm(W)
SC-NEU/176	180mm(L) x 5mm(W)
SC-NEU/177	200mm(L) x 3mm(W)
SC-NEU/178	200mm(L) x 4mm(W)
SC-NEU/179	200mm(L) x 5mm(W)
SC-NEU/180	200mm(L) x 6mm(W)
SC-NEU/181	230mm(L) x 2mm(W)
SC-NEU/182	230mm(L) x 3mm(W)
SC-NEU/183	230mm(L) x 4mm(W)
SC-NEU/184	230mm(L) x 5mm(W)

130° downward cutting

SC-NEU/198	180mm(L) x 1mm(W)
SC-NEU/199	180mm(L) x 2mm(W)
SC-NEU/200	180mm(L) x 3mm(W)
SC-NEU/201	180mm(L) x 4mm(W)
SC-NEU/202	180mm(L) x 5mm(W)
SC-NEU/203	230mm(L) x 2mm(W)
SC-NEU/204	230mm(L) x 3mm(W)
SC-NEU/205	230mm(L) x 4mm(W)
SC-NEU/206	230mm(L) x 5mm(W)

Detachable

90° upwardcutting

SC-NEU/145	180(L) x 1mm(W)
SC-NEU/146	180(L) x 2mm(W)
SC-NEU/147	180(L) x 3mm(W)
SC-NEU/148	180(L) x 4mm(W)
SC-NEU/149	180(L) x 5mm(W)
SC-NEU/150	230(L) x 2mm(W)
SC-NEU/151	230(L) x 3mm(W)
SC-NEU/152	230(L) x 4mm(W)
SC-NEU/153	230(L) x 5mm(W)

90° downward cutting

SC-NEU/163	180(L) x 1mm(W)
SC-NEU/164	180(L) x 2mm(W)
SC-NEU/165	180(L) x 3mm(W)
SC-NEU/166	180(L) x 4mm(W)
SC-NEU/167	180(L) x 5mm(W)
SC-NEU/168	230(L) x 2mm(W)
SC-NEU/169	230(L) x 3mm(W)
SC-NEU/170	230(L) x 4mm(W)
SC-NEU/171	230(L) x 5mm(W)

130° upward cutting

SC-NEU/185	180(L) x 1mm(W)
SC-NEU/186	180(L) x 2mm(W)
SC-NEU/187	180(L) x 3mm(W)
SC-NEU/188	180(L) x 4mm(W)
SC-NEU/189	180(L) x 5mm(W)
SC-NEU/190	200(L) x 3mm(W)
SC-NEU/191	200(L) x 4mm(W)
SC-NEU/192	200(L) x 5mm(W)
SC-NEU/193	200(L) x 6mm(W)
SC-NEU/194	230(L) x 2mm(W)
SC-NEU/195	230(L) x 3mm(W)
SC-NEU/196	230(L) x 4mm(W)
SC-NEU/197	230(L) x 5mm(W)

130° downward cutting

SC-NEU/207	180(L) x 1mm(W)
SC-NEU/208	180(L) x 2mm(W)
SC-NEU/209	180(L) x 3mm(W)
SC-NEU/210	180(L) x 4mm(W)
SC-NEU/211	180(L) x 5mm(W)
SC-NEU/212	230(L) x 2mm(W)
SC-NEU/213	230(L) x 3mm(W)
SC-NEU/214	230(L) x 4mm(W)
SC-NEU/215	230(L) x 5mm(W)

Adson Dissecting Forceps


Bayonet, 1-2 teeth, jaw 2.5mm at tip



SC-NEU/001c	191mm long
-------------	------------

Adson Dissecting Forceps

Bayonet, fine, serrated jaws, jaw 2.5mm at tip



SC-NEU/009	178mm long
------------	------------

Northfield Dressing Forceps

Extra light, serrated jaws



SC-NEU/014	1mm wide, 203mm long
------------	----------------------

Jefferson Dressing Forceps




SC-NEU/013	Serrated jaws, 178mm long
------------	---------------------------

Jefferson Forceps



SC-NEU/002	Fine, 1-2 teeth, 178mm long
------------	-----------------------------

Pennybacker Dissecting Forceps



SC-NEU/005	3-4 teeth, 165mm long
------------	-----------------------

Silcock Iris Forceps

102mm long



SC-NEU/006	1-2 teeth
SC-NEU/015	Serrated jaws

Cairns Dissecting Forceps

Delicate, 2-3 teeth



SC-NEU/003a	165mm long
SC-NEU/003b	216mm long

Cairns (Dott) Dressing Forceps

Serrated jaws



SC-NEU/010	Medium, 2mm wide, 216mm long
SC-NEU/012	Fine, 1mm wide, 178mm long

Dott Dissecting Forceps


Fine, 1-2 teeth



SC-NEU/004	178mm long
------------	------------

Dott Dressing Forceps


Angled



SC-NEU/018	114mm long
------------	------------

Tilley Dressing Forceps


Angled, serrated ends



SC-NEU/017	89mm long
------------	-----------

Cairns Scissors

Angled on flat



SC-NEU/020	178mm long
------------	------------

McKissock Scissors


Straight



SC-NEU/021a	146mm long
-------------	------------

Olivecrona Scissors


Angled, with all micro serrated blade



SC-NEU/024	22mm blade, 150mm shaft
------------	-------------------------

Schmieden (Taylor) Scissors


Angled, probe-pointed under blade



SC-NEU/023	146mm long
------------	------------

Sweet Scissors

Delicate



SC-NEU/025a	7mm blade, 144mm shaft
SC-NEU/030	13mm blade, 117mm shaft