


BONE FILES AND RASPS



Bone files and rasps are used for shaping, smoothing, and refining bone surfaces during orthopedic surgeries. Our instruments feature sharp, precision-cut teeth that provide efficient material removal while maintaining control and accuracy. Made from high-grade surgical steel, they offer durability, corrosion resistance, and a comfortable grip for consistent performance in demanding surgical environments.


BONE FILES

BONE FILE #10



| CODE # | DESCRIPTION |
|---------|--|
| AW69542 | Double-ended, plain serrations, 5 mm width, 7-1/8" (18 cm) |

BONE FILE #12




| CODE # | DESCRIPTION |
|---------|--|
| AW69543 | Double-ended, plain serrations, 7 and 5 mm width, 7-1/8" (18 cm) |

BONE FILE #12A




| CODE # | DESCRIPTION |
|---------|--|
| AW69544 | Double-ended, plain and cross serrations, 5 mm width, 7-1/8" (18 cm) |

BONE FILE #12CA




| CODE # | DESCRIPTION |
|---------|--|
| AW69545 | Double-ended, cross serrations, 5 and 7 mm width, 7-1/8" (18 cm) |

BONE FILE #33



| CODE # | DESCRIPTION |
|---------|---|
| AW69546 | Double-ended, plain serrations, 6 and 5 mm width 7-1/8" (18 cm) |

BONE FILE #45



| CODE # | DESCRIPTION |
|---------|--|
| AW69549 | Double-ended, plain serrations, 6 and 5 mm width, 7-1/4" (18.5 cm) |

BONE RASP #64



| CODE # | DESCRIPTION |
|---------|--|
| AW58070 | Double-ended, plain serrations, 3 and 4 mm width, 7-1/8" (18 cm) |

BONE FILE #92A


Single-ended, delicate, plain serrations.



| CODE # | DESCRIPTION |
|---------|-----------------|
| AW69530 | Angled downward |

BONE FILE #92B

Single-ended, delicate, plain serrations.



| CODE # | DESCRIPTION |
|---------|---------------|
| AW69534 | Angled upward |

NAIL RASP #93



| CODE # | DESCRIPTION |
|---------|---|
| AW69538 | Double-ended, 2 mm wide, angled, 6-3/4" (17 cm) |

BONE FILE

Single-ended, 2-sided, convex and flat blade 23 mm wide tapers to 6 mm.



| CODE # | DESCRIPTION |
|---------|-------------|
| AW46705 | 11" (28 cm) |

BONE FILE

Single-ended, bayonet shank, half round blade.



| CODE # | DESCRIPTION |
|---------|--|
| AW46650 | Fine teeth, 15 mm wide, 9-1/2" (24 cm) |

BONE RASPS

COTTLE RASP



| CODE # | DESCRIPTION |
|--------|---|
| 052130 | Straight down-cutting blade, 9 mm width, 8-1/4" (21 cm) |

AUFRICHT-GLABELLA RASP



| CODE # | DESCRIPTION |
|--------|--|
| 052165 | Up and down-cutting, 9 x 30 mm, 7-1/2" (19 cm) |


CONVERSE RASP



| CODE # | DESCRIPTION |
|--------|--|
| 052180 | Double-ended, fine and coarse serrations, 7-1/8" (18 cm) |

BONE RASP

For rapid shaping of bone.



| CODE # | DESCRIPTION |
|---------|--|
| AW69547 | Double-ended, fine and coarse cross serrations, 13 mm wide, 8-1/2" (21.5 cm) |

PUTTI BONE RASP

Double-ended. Flat blades, 18 mm wide, taper to 4 mm. One blade end is curved up and one end is curved sideways.



| CODE # | DESCRIPTION |
|---------|-------------------|
| AW46654 | 10-3/8" (26.5 cm) |

PUTTI BONE RASP

Double-ended. One round blade, 10 mm wide, tapers to 4 mm, and one half round blade, 15 mm wide, tapers to 5 mm.



| CODE # | DESCRIPTION |
|---------|---------------|
| AW46652 | 12" (30.5 cm) |

AUFRICHT RASPS



| CODE # | DESCRIPTION |
|---------|--|
| 052162 | Forward-cutting, 25 x 9 mm, 8" (20.5 cm) |
| AW58050 | Forward-cutting, 8.5 mm wide, 8-1/4" (21 cm) |
| AW58052 | Backward-cutting, 9 mm wide, 8-1/4" (21 cm) |

FOMAN RASP

Double-ended. Fine and coarse serrations. 4-sided, convex and flat blades.



| CODE # | DESCRIPTION |
|---------|------------------|
| AW46655 | 8-1/2" (21.5 cm) |

JOSEPH RASP



| CODE # | DESCRIPTION |
|--------|-----------------------------------|
| 441530 | Fine serrations, 6-1/8" (15.5 cm) |

KLEINERT-KUTZ BONE RASP

Double-ended. Fine and coarse serrations. 4-sided, convex and flat blades.



| CODE # | DESCRIPTION |
|---------|------------------|
| AW46653 | 8-7/8" (22.5 cm) |

KLEINERT-KUTZ RASP



| CODE # | DESCRIPTION |
|--------|---|
| 441555 | 4-sided tips, 7.5 mm wide, 8-5/8" (22 cm) |

LEWIS RASPS

Single-ended bayonet shank, plain serrations, 8 mm wide, 7-1/2" (19 cm).



| CODE # | DESCRIPTION |
|---------|-------------------|
| AW58072 | Fine serrations |
| AW58073 | Coarse serrations |

MALTZ RASP



| CODE # | DESCRIPTION |
|--------|---|
| 052100 | Straight down-cutting blade, 6 mm width, 7-1/8" (18 cm) |


MALTZ RASP IMPROVED



| CODE # | DESCRIPTION |
|---------|---|
| AW58062 | Straight down-cutting blade, 6 mm width, 7-1/8" (18 cm) |

POLOKOFF BONE RASPS

Double-ended, plain serrations.



| CODE # | DESCRIPTION |
|---------|---------------------------------|
| AW69540 | 6 mm wide, 7-7/8" (20 cm) |
| AW69541 | 3 and 4 mm wide, 7-7/8" (20 cm) |

RIB RASPS ALEXANDER PERIOSTEAL RASP




| CODE # | DESCRIPTION |
|---------|-------------------------------|
| AW58110 | 11 mm sharp, 8-1/8" (20.5 cm) |

FARABEUF PERIOSTEAL RASPS



| CODE # | DESCRIPTION |
|---------|-----------------------------------|
| AW58106 | Straight blade, 13 mm curved edge |
| AW58108 | Curved blade, 13 mm straight edge |

KIRMIAWION PERIOSTEAL RASP



| CODE # | DESCRIPTION |
|---------|-----------------------|
| AW58112 | Sharp, 6-3/4" (17 cm) |