




Bone tamps and impactors are designed to assist in positioning and securing bone grafts or implants during orthopedic procedures. Bone tamps provide controlled pressure to shape or compact graft material, while impactors deliver precise force to seat implants without causing damage. Manufactured from durable surgical-grade materials, these instruments offer reliability, balance, and ergonomic handling for accurate surgical execution.

BONE TAMPS AND IMPACTORS

BONE TAMPS


Cross-serrated tip, 6-1/4" (16 cm).



CODE #	DESCRIPTION
AW46893	2 mm diameter
AW46894	3 mm
AW46895	4 mm
AW46896	5 mm
AW46897	6 mm
AW46898	8 mm
AW46899	10 mm

BONE IMPACTOR

Designed to tap bone graft or bone parts into place with minimal trauma to the bone. The compactor head is nylon with a grid surface to prevent slippage.




CODE #	DESCRIPTION
AW69980	Anodized aluminum handle with replaceable nylon caps, head diameter 20 mm, 6-1/2" (16.5 cm)

Autoclavable to 250°F (121°C)

BONE CHIP PACKERS

Stainless steel handle, 9-7/8" (25 cm).



CODE #	DESCRIPTION
AW49750	3 mm x 16 mm
AW49751	2 mm x 9 mm
AW49752	6 mm x 9 mm
AW49753	6 mm x 19 mm
AW49754	1/4" (5 mm)

IMPACTORS

ENDOSCOPIC BONE TAMP

Straight, stainless steel handle, 17" (43 cm).



LUMBAR AND CERVICAL IMPACTORS AND DRIVERS



CODE #	TIP WIDTH	TIP LENGTH	SHAFT LENGTH	SHAFT DIAMETER	HANDLE LENGTH
AW55676	3/8" (10 mm)	0.16"	12.59"	0.280" (7.1 mm)	4-1/4"

CODE #	DESCRIPTION
AW55370	Impactor, standard offset, 6 mm x 12 mm, 8-1/2" (21.5 cm)
AW55371	Impactor, central canal, 7 mm x 16 mm, 8-1/2" (21.5 cm)
AW55372	Impactor, lateral narrow, 6 mm x 13 mm, 8-1/2" (21.5 cm)
AW55373	Impactor, micro offset, 4 mm x 9.5 mm, 7" (18 cm)
AW55374	Impactor, micro lateral, 4 mm x 10 mm, 7" (18 cm)
AW55375	Graft driver, lumbar, guarded, 10 mm x 13 mm, 8-1/2" (21.5 cm)
AW55376	Graft driver, anterior, guarded, 7 mm x 13 mm, 8-1/2" (21.5 cm)
AW55377	Graft driver, lumbar, unguarded, 10 mm x 13 mm, 8-1/2" (21.5 cm)
AW55378	Graft driver, anterior, unguarded, 7 mm x 13 mm, 8-1/2" (21.5 cm)