

CURETTES


Curettes are precision instruments used for scraping or debriding bone and soft tissue in orthopedic procedures. Designed with sharp, scoop-shaped tips, they allow for controlled removal of unwanted tissue while preserving surrounding structures. Manufactured from surgical-grade stainless steel, our curettes offer excellent balance, durability, and tactile feedback for precise surgical work.



BRUN CURETTES

BRUN (SPRATT) CURETTES

Oval cup, hollow handle




CODE #	TIP SHAPE	SIZE	LENGTH
AW46511	Straight	#5/0	7" (18 cm)
AW46512	Straight	#4/0	7" (18 cm)
AW46513	Straight	#3/0	7" (18 cm)
AW46515	Straight	#2/0	7" (18 cm)
AW46517	Straight	#0	7" (18 cm)
AW46519	Straight	#1	7" (18 cm)
AW46521	Straight	#2	7" (18 cm)
AW46523	Straight	#3	7" (18 cm)
AW46525	Straight	#4	7" (18 cm)
AW46529	Straight	#5	7" (18 cm)
AW46531	Straight	#6	7" (18 cm)
AW46563	Straight	#3/0	9" (23 cm)
AW46564	Straight	#2/0	9" (23 cm)
AW46565	Straight	#0	9" (23 cm)
AW46566	Straight	#1	9" (23 cm)
AW46567	Straight	#2	9" (23 cm)
AW46568	Straight	#3	9" (23 cm)
AW46569	Straight	#4	9" (23 cm)
AW46570	Straight	#5	9" (23 cm)
AW46571	Straight	#6	9" (23 cm)

CODE #	TIP SHAPE	SIZE	LENGTH
AW49110	Angled 45 degrees	#3/0	9" (23 cm)
AW49111	Angled 45 degrees	#2/0	9" (23 cm)
AW49112	Angled 45 degrees	#0	9" (23 cm)
AW49113	Angled 45 degrees	#1	9" (23 cm)
AW49114	Angled 45 degrees	#2	9" (23 cm)
AW49115	Angled 45 degrees	#3	9" (23 cm)
AW49116	Angled 45 degrees	#4	9" (23 cm)
AW49117	Angled 45 degrees	#5	9" (23 cm)
AW49118	Angled 45 degrees	#6	9" (23 cm)

BRUN (SPRATT) CURETTES


Round up, hollow handle



CODE #	TIP SHAPE	SIZE	LENGTH
AW49120	Straight	#3/0	6-3/4" (17 cm)
AW49121	Straight	#2/0	6-3/4" (17 cm)
AW49122	Straight	#0	6-3/4" (17 cm)
AW49123	Straight	#1	6-3/4" (17 cm)
AW49124	Straight	#2	6-3/4" (17 cm)
AW49125	Straight	#3	6-3/4" (17 cm)
AW49126	Straight	#4	6-3/4" (17 cm)

BRUN CURETTES


8-7/8" (22.5 cm).



CODE #	DESCRIPTION
343500	Size 000
343501	Size 00
343501	Size 0
343503	Size 1
343504	Size 2
343505	Size 3
343506	Size 4

BRUN CURETTES

Solid hexagonal handle, 9" (23 cm).



CODE #	TIP SHAPE	SIZE
AW46551	Straight	#3/0
AW46552	Straight	#2/0
AW46553	Straight	#0
AW46554	Straight	#1
AW46555	Straight	#2
AW46556	Straight	#3
AW46557	Straight	#4
AW46558	Straight	#5
AW46559	Straight	#6

CODE #	TIP SHAPE	SIZE
AW49150	Angled 45 degrees	#3/0
AW49151	Angled 45 degrees	#2/0
AW49152	Angled 45 degrees	#0
AW49153	Angled 45 degrees	#1
AW49154	Angled 45 degrees	#2
AW49155	Angled 45 degrees	#3
AW49156	Angled 45 degrees	#4
AW49157	Angled 45 degrees	#5
AW49158	Angled 45 degrees	#6

COBB CURETTES

COBB CURETTES, ALUMINUM HANDLE


All stainless steel working end and anodized aluminum, knurled handle, 11" (28 cm).



CODE #	TIP SHAPE	SIZE	CATALOG #	TIP SHAPE	SIZE
AW46433	Angled	#3/0	AW46422	Straight	#3/0
AW46434	Angled	#2/0	AW46423	Straight	#2/0
AW46435	Angled	#0	AW46424	Straight	#0
AW46436	Angled	#1	AW46425	Straight	#1
AW46437	Angled	#2	AW46426	Straight	#2
AW46438	Angled	#3	AW46427	Straight	#3
AW46439	Angled	#4	AW46428	Straight	#4
AW46440	Angled	#5	AW46429	Straight	#5
AW46441	Angled	#6	AW46430	Straight	#6

COBB CURETTES, STAINLEAW STEEL HANDLE

All stainless steel working end and stainless steel, knurled handle, 11" (28 cm).




CODE #	TIP SHAPE	SIZE	CODE #	TIP SHAPE	SIZE
AW46433S	Angled	#3/0	AW46422S	Straight	#3/0
AW46434S	Angled	#2/0	AW46423S	Straight	#2/0
AW46435S	Angled	#0	AW46424S	Straight	#0
AW46436S	Angled	#1	AW46425S	Straight	#1
AW46437S	Angled	#2	AW46426S	Straight	#2
AW46438S	Angled	#3	AW46427S	Straight	#3
AW46439S	Angled	#4	AW46428S	Straight	#4
AW46440S	Angled	#5	AW46429S	Straight	#5
AW46441S	Angled	#6	AW46430S	Straight	#6

SPECIALTY BONE CURETTES


EPSTEIN CURETTES

Reverse angle, hollow handle, 8" (20.5 cm).



CODE #	DESCRIPTION
AW49580	#000
AW49582	#00
AW49584	#0

SCOVILLE CURETTES



CODE #	DESCRIPTION
AW49574	Straight, 10" (25.5 cm)
AW49576	Down angle, 10" (25.5 cm)
AW49578	Up angle, 10" (25.5 cm)

JANSEN CURETTE



CODE #	DESCRIPTION
AW29632	Double-ended, oval cups, 6" (15 cm)

MARTINI BONE CURETTE



CODE #	DESCRIPTION
AW29620	Double-ended, 4 mm and 5 mm round cups, 5-1/2" (14 cm)

HATFIELD BONE CURETTES

Stainless steel with satin finish. Oval shaped tips, overall length 9" (23 cm). Complete set of nine curettes



CODE #	TIP SHAPE	SIZE	TIP LENGTH
584160	000	2.5 mm	4.0 mm
	Straight	2.8 mm	4.6 mm
	Straight	3.3 mm	5.0 mm
	Straight	3.8 mm	6.0 mm
	Straight	4.8 mm	7.6 mm
	Straight	5.8 mm	9.2 mm
	Straight	7.2 mm	11.6 mm
	Straight	8.7 mm	12.8 mm
	Straight	10.5 mm	15.5 mm

VOLAWANN CURETTES

Double-ended, oval and round cups.



CODE #	DESCRIPTION
AW46560	5" (12.5 cm)
AW46561	8" (20.5 cm)
AW46562	6-1/2" (16.5 cm)


WILLIGER BONE CURETTE



CODE #	DESCRIPTION
AW29630	Double-ended, 3 mm and 4 mm oval cups, 5-1/2" (14 cm)

VOLAWANN CURETTES


Double-ended, oval cups.



CODE #	DESCRIPTION
AW49200	Double-ended, straight, 5" (12.5 cm)
AW49202	Double-ended, straight, 5 1/2" (14 cm)
AW46550	Double-ended, oval cups, 6 1/2" (16.5 cm)
AW49204	Double-ended, straight, 8" (20.5 cm)

RING CURETTES


Hollow handle with fenestrated cup.



CODE #	DESCRIPTION
AW55680	Size 0, small, 10" (25.5 cm)
AW55684	Size 0, small, 11-1/2" (29 cm)
AW55681	Size 2, medium, 10" (25.5 cm)
AW55685	Size 2, medium, 11-1/2" (29 cm)
AW55682	Size 6, large, 10-1/4" (26 cm)
AW55686	Size 6, large, 12" (30.5 cm)

CONE RING CURETTES


Knurled handle and fenestrated cup, 9" (23 cm).



CODE #	TIP SHAPE	SIZE	
AW49450	Straight	#1	8 mm diameter
AW49452	Straight	#2	6 mm diameter
AW49454	Straight	#3	3 mm diameter
AW49456	Angled 35 degrees	#1	8 mm diameter
AW49458	Angled 35 degrees	#2	6 mm diameter
AW49460	Angled 35 degrees	#3	3 mm diameter

RAY PITUITARY CURETTES

7-1/2" (19 cm).




CODE #	TIP SHAPE	SIZE
AW49403	Vertical angle 45 degrees	3 mm
AW49404	Vertical angle 45 degrees	4 mm
AW49405	Vertical angle 45 degrees	5 mm
AW49413	Vertical angle 90 degrees	3 mm
AW49414	Vertical angle 90 degrees	4 mm
AW49415	Vertical angle 90 degrees	5 mm
AW49423	Horizontal angle 45 degrees	3 mm
AW49424	Horizontal angle 45 degrees	4 mm
AW49425	Horizontal angle 45 degrees	5 mm
AW49433	Horizontal angle 90 degrees	3 mm
AW49434	Horizontal angle 90 degrees	4 mm
AW49435	Horizontal angle 90 degrees	5 mm

ENDOSCOPIC CURETTES



ENDO BRUN CURETTES

Endoscopic with phenolic handle



CODE #	TIP WIDTH	TIP	TIP LENGTH	SHAFT LENGTH	SHAFT DIAMETER	HANDLE LENGTH	OVERALL LENGTH
AW55532	Angled 40°	2 x 4 mm cup	2.7"	9.3"	0.285"	5"	17" (43 cm)
AW55523	Straight	3 x 6 mm cup	2.7"	9.3"	0.285"	5"	17" (43 cm)
AW55533	Angled 40°	3 x 6 mm cup	2.7"	9.3"	0.285"	5"	17" (43 cm)
AW55543	Angled 40°	3 x 6 mm cup	3-1/2"	10"	0.45"	6-1/2"	20" (51 cm)
AW55556	Angled 30°	6 x 9 mm cup	3-1/2"	10"	0.45"	6-1/2"	20" (51 cm)


ENDO EPSTEIN CURETTES



CODE #	TIP WIDTH	TIP	TIP LENGTH	SHAFT LENGTH	SHAFT DIAMETER	HANDLE LENGTH	OVERALL LENGTH
AW55593	Reverseangled 40°	3 x 6 mm	2.7"	9.3"	0.285"	5"	17" (43 cm)
AW55586	Reverseangled 30°	6 x 9 mm	3-1/2"	10"	0.45"	6-1/2"	20" (51 cm)

ENDO CONE RING CURETTES


Endoscopic with phenolic handle



CODE #	TIP WIDTH	TIP	TIP LENGTH	SHAFT LENGTH	SHAFT DIAMETER	HANDLE LENGTH	OVERALL LENGTH
AW55573	Angled 40°	3 mm	2.7"	9.3"	0.285"	5"	17" (43 cm)
AW55566	Angled 30°	6 mm	10"	10"	0.45"	6-1/2"	20" (51 cm)

ENDO SCOVILLE CURETTE

Endoscopic with phenolic handle.



CODE #	TIP WIDTH	TIP	TIP LENGTH	SHAFT LENGTH	SHAFT DIAMETER	HANDLE LENGTH	OVERALL LENGTH
AW55510	Straight	9 x 15 mm	3-1/2"	10"	0.45"	6-1/2"	20" (51 cm)

SMARTRACK CURETTE TRAYS

BRUN SMARTRACK FOR 7" CURETTES

Aluminum and stainless steel.



CODE #	DESCRIPTION
AW98160	Holds 12 7" AWedic® Brun Currettes

COBB SMARTRACK FOR 11" INSTRUMENTS

Aluminum and stainless steel.



CODE #	DESCRIPTION
AW98170	Holds nine 11" AWedic® Cobb Currettes, Elevators or Gouges

HIBBS AND BRUN SMARTRACK FOR 9" INSTRUMENTS

Aluminum and stainless steel.



CODE #	DESCRIPTION
AW98180	Holds 12 9" AWedic® Brun Currettes, or Hibbs Chisels, Gouges or Osteotomes