

MALLETS

Mallets are used to deliver controlled force to chisels, osteotomes, and other orthopedic instruments during surgical procedures. Our mallets are designed for balance, precision, and comfort, ensuring effective impact with minimal vibration. Constructed from high-quality surgical materials, they are available in various weights and head designs to suit different surgical needs.



LIGHT & MEDIUM WEIGHT MALLET

NARROW TIP MALLET

Useful for hand surgery and in areas where access is limited.



Autoclavable to 290°F (143°C)

CODE #	WEIGHT	HEAD WEIGHT	HEAD DIAMETER	OVERALL LENGTH
AW49764	7 oz (198 g)	3 oz (85 g)	18 x 18 mm (square contact surface) 7 mm (round contact surface)	8" (20.5 cm)

MINI MALLET

Solid stainless steel.



CODE #	WEIGHT	HEAD WEIGHT	HEAD DIAMETER	OVERALL LENGTH
AW49765	4 oz (113 g)	3 oz (85 g)	20 mm	7" (18 cm)

SOLID COPPER HEAD MALLET

Solid stainless steel.



CODE #	WEIGHT	HEAD DIAMETER	OVERALL LENGTH
051510	4 oz (113 g)	3/4" (20 mm)	6-1/2" (16.5 cm)

LUCAE BONE MALLET

Solid stainless steel handles.



CODE #	WEIGHT	HEAD	HEAD WEIGHT	HEAD DIAMETER	OVERALL LENGTH
AW49761	7 oz (198 g)	Stainless steel	3 oz (85 g)	19 mm	7-1/2" (19 cm)
AW46660	9 oz (260 g)	Brass	5 oz (142 g)	25 mm	8" (20.5 cm)

MALLET #49

Lead filled and faced



CODE #	WEIGHT	HEAD WEIGHT	HEAD DIAMETER	OVERALL LENGTH
AW46656	7 oz (198 g)	5 oz (142 g)	22 mm	7-3/4" (19.5 cm)

MALLET CAPS



CODE #	DESCRIPTION
AW46662P	Pair of nylon caps for AW46662, 1" (2.5 cm)

NYLON HEAD MALLETS



CODE #	WEIGHT	HEAD DIAMETER	OVERALL LENGTH
051520	8 oz (227 g)	1" (25 mm)	7-1/2" (19 cm)
051525	11 oz (312 g)	1-1/4" (32 mm)	7-5/8" (19.5 cm)

MALLET

Lightweight. Includes two replaceable nylon caps, anodized aluminum handle.




CODE #	WEIGHT	HEAD WEIGHT	HEAD DIAMETER	OVERALL LENGTH
AW46662	7 oz (198 g)	4 oz (113 g)	25 mm	7-1/2" (19 cm)

NYLON HEAD MALLET FACES



CODE #	DESCRIPTION
051521	Extra 1" nylon face for 051520
051526	Extra 1-1/2" nylon face for 051525

PECK MALLET



CODE #	WEIGHT	OVERALL LENGTH
441420	7.7 oz (218 g)	7-7/8" (20 cm)

COLLIN MALLET

Solid stainless steel



CODE #	WEIGHT	HEAD WEIGHT	HEAD DIAMETER	OVERALL LENGTH
AW49760	7 oz (198 g)	5 oz (142 g)	30 mm	8" (20.5 cm)


MEAD MALLET



CODE #	WEIGHT	HEAD WEIGHT	HEAD DIAMETER	OVERALL LENGTH
AW46657	11 oz (315 g)	6 oz (170 g)	20 mm	7-1/2" (19 cm)

COMBINATION MALLET

With one nylon cap and one solid stainless steel side.



CODE #	WEIGHT	HEAD WEIGHT	HEAD DIAMETER	OVERALL LENGTH
AW46686	9 oz (260 g)	5 oz (142 g)	25 mm	7-1/2" (19 cm)
AW46687	12 oz (340 g)	8 oz (227 g)	32 mm	7-1/2" (19 cm)

GERZOG BONE MALLET


Lead filled and faced



CODE #	WEIGHT	HEAD WEIGHT	HEAD DIAMETER	OVERALL LENGTH
AW46661	10 oz (285 g)	7 oz (200 g)	25 mm	7-1/2" (19 cm)

BERGMAN MALLET

Aluminum alloy handle.



CODE #	WEIGHT	HEAD WEIGHT	HEAD DIAMETER	OVERALL LENGTH
AW49766	15 oz (425 g)	9 oz (255 g)	30 mm	8-7/8" (22.5 cm)
AW49767	1 lb 1 oz (482 g)	11 oz (312 g)	45 mm	9-1/2" (24 cm)

CLOWARD-STYLE MALLET

Solid stainless steel.



CODE #	WEIGHT	HEAD WEIGHT	HEAD DIAMETER	BALL PEEN DIAMETER	OVERALL LENGTH
AW55400	8 oz (227 g)	4 oz (113 g)	20.5 mm	17.5 mm	7-1/2" (19 cm)
AW55401	16 oz (454 g)	8 oz (227 g)	25 mm	22 mm	9" (23 cm)

COMBINATION MALLET CAPS



CODE #	DESCRIPTION
AW46686P	Nylon cap for AW46686, 1" (2.5 cm)
AW46687P	Nylon cap for AW46687, 1-1/4" (3.2 cm)

ORTHOPEDIC MALLET

Solid stainless steel.




CODE #	HANDLE LENGTH	WEIGHT	HEAD WEIGHT	HEAD DIAMETER	OVERALL LENGTH
AW49763	Short	1 lb 2 oz (510 g)	8 oz (227 g)	30.5 mm	7" (18 cm)
AW46666	Normal	2 lbs (900 g)	1 lb (455 g)	32 mm	10-3/4" (27.5 cm)

HEAVYWEIGHT MALLET

KIRK MALLET

Brass filled with stainless steel face.



Head is brass filled and should only be used to strike flat surfaces

CODE #	WEIGHT	HEAD WEIGHT	HEAD DIAMETER	OVERALL LENGTH
AW46658	1 lb 7 oz (650 g)	1 lb 4 oz (570 g)	38 mm	8" (20.5 cm)

ORTHOPEDIC MALLET

Solid stainless steel.



CODE #	HANDLE LENGTH	WEIGHT	HEAD WEIGHT	HEAD DIAMETER	OVERALL LENGTH
AW49762	Short	2 lbs (900 g)	1 lb 5 oz (595 g)	38 mm	7-1/2" (19 cm)
AW46667	Normal	3 lbs (1360 g)	2 lbs (900 g)	50 mm	11" (28 cm)

HEATH MALLET


Heavyweight, solid stainless steel.



CODE #	WEIGHT	HEAD WEIGHT	HEAD DIAMETER	OVERALL LENGTH
AW46663	2 lbs (907 g)	1 lb 4 oz (570 g)	40 mm	7-1/2" (19 cm)

PHENOLIC HANDLE MALLET

Grooved phenolic handle for comfortable and safe grip



Autoclavable to 290°F (143°C)

CODE #	WEIGHT	HEAD WEIGHT	HEAD DIAMETER	OVERALL LENGTH
AW46688	1 lb 11 oz (765 g)	(765 g) 1 lb (455 g)	32 mm	11" (28 cm)
AW46689	2 lbs 9 oz (1162 g)	2 lbs (900 g)	45 mm	11" (28 cm)

REPERCUSSION FREE MALLET

Also known as a dead blow mallet. One nylon cap and one stainless steel side. Grooved phenolic handle for comfortable grip



Autoclavable to 250°F (121°C)

CODE #	WEIGHT	HEAD WEIGHT	HEAD DIAMETER	OVERALL LENGTH
AW46690	1 lb 5 oz (600 g)	10 oz (285 g)	30 mm	11" (28 cm)

REPERCUSSION FREE MALLET CAP



CODE #	DESCRIPTION
AW46690P	Cap only for AW46690

BONE CUTTING FORCEPS


CLEVELAND BONE CUTTING FORCEPS

For small bones.



CODE #	DESCRIPTION
AW47270	5-3/4" (15 cm)
AW47272	6-1/2" (17 cm)

MCINDOE BONE CUTTING FORCEPS




CODE #	DESCRIPTION
AW47546	Angled, very delicate, 7" (18 cm)

LISTON BONE CUTTING FORCEPS




CODE #	DESCRIPTION
AW47398	Straight, 5-1/2" (14 cm)
AW47410	Angled on flat, 5-1/2" (14 cm)
AW47400	Straight, 6-3/4" (17 cm)
AW47404	Angled on flat, 6-3/4" (17 cm)
AW47402	Straight, 7-7/8" (20 cm)
AW47412	Angled, 7-7/8" (20 cm)
AW47406	Straight, 8-5/8" (22 cm)
AW47408	Angled on flat, 8-1/2" (22 cm)



RUSKIN-LISTON BONE CUTTING FORCEPS

Delicate blades.



CODE #	DESCRIPTION
AW47555	Straight, 5-5/8" (15 cm)
AW47557	Angled, 5-1/2" (14 cm)
AW47540	Straight, 7-1/8" (18 cm)
AW47541	Angled on flat, 7" (18 cm)


STILLE-LISTON BONE CUTTING FORCEPS

10-1/2" (27 cm).




CODE #	DESCRIPTION
AW47366	Straight
AW47370	Angled on flat

STILLE-HORSLEY BONE CUTTING FORCEPS



CODE #	DESCRIPTION
AW47374	Angled, 9-3/8" (24 cm)

RUSKIN-LISTON BONE CUTTING FORCEPS



CODE #	DESCRIPTION
AW47550	Straight, 6" (15 cm)
AW47538	Straight, 7-1/4" (19 cm)
AW47542	Angled, 7-1/4" (19 cm)