

OSTEOTOMES


Osteotomes are used to cut or reshape bone with controlled force during orthopedic surgeries. Designed with sharp, tapered blades, they allow for precise bone sectioning while minimizing trauma to surrounding tissue. Made from high-strength surgical steel, our osteotomes offer excellent durability, balance, and grip for safe and effective use in a variety of procedures.



LAMBOTTE


MINI-LAMBOTTE OSTEOTOMES, STRAIGHT 5"

Straight, without calibration lines, 5" (12.5 cm).



CODE #	DESCRIPTION
AW46336	1/8" (4 mm)
AW46338	1/4" (6 mm)
AW46340	5/16" (8 mm)
AW46342	3/8" (10 mm)
AW46344	1/2" (12 mm)
AW46343	9/16" (14 mm)
AW46345	5/8" (15 mm)
AW46346	11/16" (18 mm)
AW46347	3/4" (19 mm)
AW46335	13/16" (20 mm)

LAMBOTTE OSTEOTOMES, STRAIGHT 7"



CODE #	DESCRIPTION
AW46348	1/8" (4 mm)
AW46349	3/16" (5 mm)
AW46350	1/4" (6 mm)
AW46358	5/16" (8 mm)
AW46351	3/8" (10 mm)
AW46352	1/2" (12 mm)
AW46353	5/8" (16 mm)
AW46354	3/4" (19 mm)
AW46355	7/8" (22 mm)
AW46356	1" (25 mm)

LAMBOTTE-STYLE OSTEOTOMES


flat thin blade, 5" (12.5 cm)



CODE #	DESCRIPTION
051082	2 mm
051084	4 mm
051086	6 mm
051088	8 mm
051090	10 mm
051092	12 mm
051095	15 mm


LAMBOTTE OSTEOTOMES, STRAIGHT 9"

Straight, without calibration lines, 9" (23 cm).



CODE #	DESCRIPTION
AW46304	1/4" (6 mm)
AW46306	1/2" (13 mm)
AW46310	3/4" (19 mm)
AW46314	1" (25 mm)
AW46318	1-1/4" (32 mm)
AW46322	1-1/2" (38 mm)


LONG BEVEL OSTEOTOMES



CODE #	DESCRIPTION
AW69940	Straight, 10 mm
AW69946	Straight, 20 mm
AW69950	Straight, 25 mm
AW69970	Curved, 5 mm
AW69975	Curved, 10 mm

LAMBOTTE OSTEOTOMES, CURVED 9"

Curved, without calibration lines, 9" (23 cm).



CODE #	DESCRIPTION
AW46328	1/4" (6 mm)
AW46329	1/2" (12 mm)
AW46330	3/4" (19 mm)
AW46331	1" (25 mm)
AW46332	1-1/4" (32 mm)
AW46333	1-1/2" (38 mm)

SPECIALTY OSTEOTOMES

Flat thin blade, 6" (15 cm).



CODE #	DESCRIPTION
051102	2 mm
051103	3 mm
051104	4 mm
051105	5 mm
051106	6 mm
051107	7 mm
051110	10 mm

MINI OSTEOTOMES


Curved, 5" (12.5 cm).



CODE #	DESCRIPTION
441207	4 mm
441208	6 mm

HOKE HOKE OSTEOTOMES

5-1/4" (13.5 cm).




CODE #	TIP	TIP WIDTH
AW46497	Straight	1/8" (3 mm)
AW46498	Straight	3/16" (4 mm)
AW46499	Straight	1/4" (6 mm)
AW46500	Straight	5/16" (8 mm)
AW46501	Straight	3/8" (10 mm)
AW46502	Straight	1/2" (13 mm)
AW46503	Straight	5/8" (16 mm)
AW46504	Straight	3/4" (19 mm)

CODE #	TIP	TIP WIDTH
AW46700	Curved	1/8" (3 mm)
AW46701	Curved	3/16" (4 mm)
AW46702	Curved	1/4" (6 mm)
AW46703	Curved	5/16" (8 mm)
AW46704	Curved	3/8" (10 mm)
AW46706	Curved	1/2" (13 mm)
AW46708	Curved	5/8" (16 mm)
AW46709	Curved	3/4" (19 mm)

SILVER

SILVER OSTEOTOMES

Single guard, 7" (18 cm)




CODE #	DESCRIPTION
AW58100	Straight
AW58102	Guard, right
AW58104	Guard, left

SHEEHAN

SHEEHAN OSTEOTOMES

Straight, 6-1/4" (16 cm)



CODE #	DESCRIPTION
AW46268	2 mm
AW46270	4 mm
AW46272	6 mm
AW46274	8 mm
AW46276	10 mm
AW46278	12 mm
AW46280	15 mm

U.S. Army

U.S. ARMY PATTERN OSTEOTOMES

7" (18 cm)




CODE #	DESCRIPTION
AW46284	6 mm
AW46286	12 mm
AW46288	18 mm

STILLE

STILLE OSTEOTOMES


Straight, 8" (20.5 cm).



CODE #	DESCRIPTION
AW49002	10 mm
AW49004	14 mm
AW49006	18 mm
AW49008	22 mm

HIBBS


HIBBS OSTEOTOMES



CODE #	DESCRIPTION
AW1011	

HIBBS ENDOSCOPIC OSTEOTOME

Straight tip, 18-1/2" (47 cm).




CODE #	TIP WIDTH	TIP LENGTH	SHAFT LENGTH	SHAFT DIAMETER	HANDLE LENGTH	HANDLE LENGTH
AW55642	3/8" (10 mm)	3-1/2"	10"	0.43" (11 mm)	10-3/4"	4-1/4"

SMITH PETERSON

SMITH PETERSON OSTEOTOME


8" (20.5 cm).



CODE #	tip	TIP WIDTH
AW46448	Straight	1/4"
AW46449	Straight	3/8"
AW46450	Straight	1/2"
AW46451	Straight	5/8"
AW46452	Straight	3/4"
AW46453	Straight	1"
AW46454	Straight	1-1/4"
AW46457	Curved	1/4"
AW46458	Curved	3/8"
AW46459	Curved	1/2"
AW46460	Curved	5/8"
AW46461	Curved	3/4"
AW46462	Curved	1"
AW46463	Curved	1-1/4"

INTERCHANGEABLE OSTEOTOME, CHISEL AND GOUGE SET

INTERCHANGEABLE OSTEOTOME, CHISEL AND GOUGE SET




CODE #	DESCRIPTION
AW49630	Complete set with handle, key and all blades
AW49632	Phenolic handpiece only, 7-3/8" (18.5 cm)
AW49634	Key only
AW49668	Osteotome blade, 5 mm
AW49672	Osteotome blade, 10 mm
AW49676	Osteotome blade, 16 mm
AW49680	Osteotome blade, 25 mm
AW49652	Chisel blade, 5 mm
AW49656	Chisel blade, 10 mm
AW49640	Gouge blade, 1" wide, 50 mm radius
AW49644	Gouge blade, 1" wide, 60 mm radius
AW49648	Gouge blade, 1" wide, 70 mm radius

COBB

COBB OSTEOTOME

All stainless steel, 11" (28 cm).



CODE #	tip	TIP WIDTH
AW49350	Straight	1/4"
AW49351	Straight	3/8"
AW49352	Straight	1/2"
AW49353	Straight	5/8"
AW49354	Straight	3/4"
AW49355	Straight	1"
AW49356	Straight	1-1/4"
AW49357	Curved	1/4"
AW49358	Curved	3/8"
AW49359	Curved	1/2"
AW49360	Curved	5/8"
AW49361	Curved	3/4"
AW49362	Curved	1"
AW49363	Curved	1-1/4"