



Our orthopedic screwdrivers are designed for precise control and secure handling of all standard implant screws. With ergonomic handles and torque-optimized shafts, they allow for accurate insertion and removal without stripping or slippage. Available in a variety of tip designs and sizes, each screwdriver ensures compatibility with a wide range of orthopedic systems.


SCREWDRIVERS

HEX HEAD HEXAGONAL SCREWDRIVERS




CODE #	DESCRIPTION
AW47952	2.5 mm screwdriver with sleeve, 7-7/8" (20 cm)
AW48368	2.5 mm screwdriver with notch, 8-1/4" (21 cm)
AW48349	3.5 mm screwdriver with notch, 9-5/8" (24.5 cm)
AW47955	Screw holding sleeve only for AW47952, 2" (5 cm)
AW48365	Screw holding sleeve only for AW48368, 3" (8 cm)
AW47983	Screw holding sleeve only for AW48349, 5" (12.5 cm)
AW47956	Screwdriver only for AW47952

LARGE HEXAGONAL SCREWDRIVER (PELVIC)



CODE #	DESCRIPTION
AW48353	3.5 mm screwdriver with notch, 11-3/4" (30 cm)
AW48364	Screwdriver only, 2.5 mm with notch, 11" (27 cm)
AW47983	Sleeve only for AW48353, 5" (12.5 cm)

HEXAGONAL SCREWDRIVER SHAFT (PELVIC)




CODE #	DESCRIPTION
AW47951	2.5 mm with notch, quick coupling end, 6-1/2" (16.5 cm)
AW47953	2.5 mm without notch, quick coupling end, 4" (10 cm)
AW47961	3.5 mm with notch, quick coupling end, 6-1/2" (16.5 cm)
AW48365	Sleeve only for AW47951, 3" (8 cm)
AW47983	Sleeve only for AW47961, 5" (12.5 cm)

SCREWDRIVER SLEEVES



CODE #	DESCRIPTION
AW48365	Sleeve only for AW47951, AW48348, AW48364, AW48368, 3" (8 cm)
AW47983	Sleeve only for AW48349, AW48353 and AW47961, 5" (12.5 cm)

LARGE HEXAGONAL SCREWDRIVER



CODE #	DESCRIPTION
AW48252 3	3.5 mm, quick coupling end, 3-7/8" (10 cm)

LARGE HEX SCREWDRIVER WITH UNIVERSAL JOINT



CODE #	DESCRIPTION
AW48396	3.5 mm tip, 11-5/8" (29.5 cm)

HEXAGONAL SCREWDRIVER SHAFT



CODE #	DESCRIPTION
AW48348	2.5 mm with notch, quick coupling end, 5-1/2" (14 cm)
AW48365	Sleeve only

LARGE HEXAGONAL SCREWDRIVER



CODE #	DESCRIPTION
AW47954	3.5 mm, quick coupling end, 3-7/8" (10 cm)

SCREWDRIVER




CODE #	DESCRIPTION
AW48670	1.5 mm hexagonal, 7-7/16" (200 mm)
AW48350	2.5 mm hexagonal, 9-27/32" (250 mm)
AW48352	3.5 mm hexagonal, 9-27/32" (250 mm)
AW48354	Single slot, 9-27/32" (250 mm)
AW48356	Cruciform, 9-27/32" (250 mm)
AW48358	Phillips, 9-27/32" (250 mm)

MINI HEXAGONAL SCREWDRIVER



CODE #	DESCRIPTION
AW47950	1.5 mm, mini quick coupling end, sleeve included, 2-1/4" (5.5 cm)
AW47950B	Shaft only, 1.5 mm with notch, mini quick coupling end, 2" (5 cm)
AW47959	Screw holding sleeve only

SINGLE SLOT SCREWDRIVER



CODE #	DESCRIPTION
AW48355	Phenolic handle, 6-3/4" (17 cm)

CRUCIFORM HEAD CRUCIFORM SCREWDRIVER



CODE #	DESCRIPTION
AW47957	1.5 / 2.0 mm screwdriver with sleeve, 7-1/2" (19 cm)
AW47959	Screw holding sleeve only for AW47957
AW47958	Screwdriver only for AW47957

MINI CRUCIFORM SCREWDRIVER BIT



CODE #	DESCRIPTION
AW45004	1.5 / 2 mm, mini quick coupling end, 2-1/2" (6.5 cm)

T-HANDLE LARGE HEX SCREWDRIVER W/ T-HANDLE

CODE #	DESCRIPTION
AW48395	3.5 mm, quick coupling end, 6-7/8" (17.5 cm)

DEPTH GAUGE SCREW DEPTH GAUGES

Used to determine screw length required during internal fracture fixation (mm graduated). Snaps apart for easy cleaning. IMPORTANT NOTE: When not in use, retract measuring end inside gauge for probe protection



CODE #	DESCRIPTION
AW47934	Mini, up to 30 mm, 4-3/4" (12 cm)
AW47936	Small, up to 50 mm, 6-1/4" (16 cm)
AW47940	Large, up to 100 mm, 8-3/4" (22 cm)
AW47935	Extra-long, up to 120 mm, 11" (28 cm)

INTERCHANGEABLE SCREWDRIVER SET

Screwdriver Sets



CODE #	DESCRIPTION
AW48473	Screwdriver handle. Allows K-Wire and pins to be inserted with Jacobs Universal Chuck AW48478. Cannulation 1/4" (6.3 mm), includes key, 5-1/2" (14 cm)
AW48474	Single slot, 5-1/2" (14 cm)
AW48475	Cruciform, 5-1/2" (14 cm)
AW48476	Phillips, 5-1/2" (14 cm)
AW48477	Small Phillips, 5-1/2" (14 cm)
AW48478	Jacobs Universal Chuck. Allows K-Wire and pins to be inserted with Jacobs Chuck. Cannulation 1/4" (6.3 mm), includes key, 3" (7.5 cm)
AW48478P	Key for AW48478 Jacobs Universal Chuck
AW48479	3.5 mm diameter, 5-1/2" (14 cm)
AW48480	2.5 mm diameter, 5-1/2" (14 cm)

SCREW HOLDING FORCEPS

SELF-RETAINING SCREW HOLDING FORCEPS



CODE #	DESCRIPTION
AW47930	3-1/2" (9 cm)

SELF-RETAINING SCREW HOLDING FORCEPS



CODE #	DESCRIPTION
AW45208	For 1.5 mm to 2.7 mm screws
AW45210	For 3.5 mm to 6.5 mm screws